

Aim of the Experiment:

To find the diameter and determine the cross sectional area of a given wire by using screw gauge.

Apparatus Required:

- Wire
- Screw gauge

Theory:

Pitch of a screw is defined as the distance between the two of its consecutive threads as measured along the axis of the screw.

$$\text{Pitch} = \frac{\text{distance moved by the screw}}{\text{no. of rotations given}}$$

$$= \frac{5 \text{ mm}}{5} = 1 \text{ mm}$$

Least Count: It is the least count / measurement which can be performed with the help of screw.

$$\text{Least count} = \frac{\text{Pitch}}{\text{Total no. of division of circumference}}$$

$$= \frac{1 \text{ mm}}{100}$$

Teacher's Signature _____

Expt. No. 01

$$= \frac{0.64 \text{ mm}}{100} \\ = 0.001 \text{ cm}$$

Area of Cross section = πr^2

Zero error :- If the linear scale of '0' does not coincide with the circular scale of '0' the zero error obtained

$$\text{Zero error} = \pm a \times \text{L.C.} =$$

→ Negative zero error =

Procedure :-

- (i) The pitch and the least count of the screw gauge are determined by the method of given in the theory.
- (ii) The wire is placed in the gap and the screw head is turned till the wire is touched by the two jaws. In their position in the main scale two jaws. In their position in the main scale reading to be noted.
- (iii) The Circular scale is noted down where the reference line coincide with the circular scale.
- (iv) The main scale reading and circular scale reading are found out and their sum gives the diameter of the wire.

Teacher's Signature _____

Expt. No. 01

(v) The above procedure is repeated for 10 times and their mean value is calculated to get diameter of the wire.

(w) The radius can be calculated by the formula $r = d/2$

Observation: To find the diameter of the given wire (d)

Result:- The Area of cross-section of the given wire is $0.2,098 \text{ cm}^2$

Name - Sujoyo R. Mohanta


Roll.no. - F211 5901 0060

Branch - Mining Engg.

Semester - 2nd Sem

Submission Date - 29/03/2022

Teacher's Signature


29/03/2022

Aim of the Experiment:

To find the volume of a Rs 2 coin using screw gauge.

Apparatus Required:

→ Screw gauge

→ 20 rupees coin

Theory:

Thickness Pitch scale reading + Circular scale reading
 Pitch scale reading (P.S.R) = Pitch (P) × No. of Complete
 rotations (N)

$$\therefore P.S.R = P \times N$$

Circular scale reading (C.S.R) = Difference between
 initial circular scale coincidence (I) and Final circular
 scale coincidence (F) × Least count

$$\therefore C.S.R = (I - F) \times L.C$$

$$\text{Hence Diameter} = [P \times N + (I - F) \times L.C]$$

If S be the area of the Rs 2 coin and t be its thickness then the volume ' V ' of irregular Rs 2 coin is given by $V = St \text{ cm}^3$.

Area of Rs-2 coin from graph paper

[{ No. of big square divisions (N) × Area of each big square division (A) } + { No. of small square division (n) × Area of each small square division (a) }]

Teacher's Signature _____

Qn.

$S = (NA + no) \text{ cm}^2$
 Thickness of R2 coin = Pitch scale reading + Circular scale reading

Qn.

$Z = P \times N + (I - F) \times L.C$
 where P = Pitch of the screw gauge, N = no. of complete rotation, I = Initial circular scale coincidence, F = Final coincidence, L.C = least count.

Determination of pitch and L.C

10 division of linear scale = 1cm / 0.5cm

1 division of linear scale = 0.1cm / 0.05cm

Here, 1 complete rotation = 0.1cm / 0.05cm

So, pitch = 0.1cm or 0.05cm

Least count = $\frac{\text{Pitch}}{\text{Total no. of divisions in the circular scale}}$

$$L.C = \frac{0.1 \text{ cm}}{100} = 0.001 \text{ cm or } \frac{0.05 \text{ cm}}{100} = 0.0005 \text{ cm}$$

Procedure:-

- The given irregular R2 coin is placed at ten different places on the graph paper and its outline is drawn at each location with the help of a pencil.

Teacher's Signature _____

Exp. No. 07

- (vi) The outlines are numbered as figure 1, figure 2 etc.
- (vii) The total no. of big square divisions and small square divisions enclosed inside the outline boundary are counted.
- (viii) The small square along the boundary is counted as one if the line covers half or more than half of its area.
- (ix) All the readings are noted in the table and the area of each figure is found out. The mean area is also counted and recorded.
- (x) The pitch and L.C. of the screw gauge are determined and by placing the Rs 2 coin in between the two jaws F and N are found out.
- (xi) All these readings are recorded in a tabular form and the thickness of the Rs 2 coin is determined. This procedure is repeated for at least 30 observations and the mean thickness is noted.
- (xii) By using the formula given in the theory, the volume of Rs 2 coin is calculated.

Observation:

Zero error is negative zero error = _____

Tabulation-I
(Thickness of Rs 2 coin) (t)

Sl No.	M.S.R (t)	Head Coincidence		C.S.R. <small>(in mm (n x LC))</small>	Total (M.S.R + C.S.R)	Mean
	in mm	Observed	Corrected			
1.	1	25	10.5	0.25	1.25	
2	1	15	9.5	0.15	1.15	

Teacher's Signature _____

Result: /

The volume of Rs 2 coin is found

AL

Name - Sujoyo Ranjan Mahanta


Regd no. - F21159010060

Branch - Mining Engg.

Semester - 2nd Sem.

Submission Date - 07/04/2022

Teacher's Signature


07/04/22

Aim of the Experiment ~ 3

To find the volume of a solid cylinder by using a vernier calliper.

Apparatus Required :-

Slide Calliper, Solid Cylinder.

Theory :-

Length of the cylinder = Main scale reading + vernier scale reading.

Vernier scale reading = Vernier coincidence \times least count

Tabulation of least count:

$$\begin{aligned} \text{Least Count} &= \frac{1 \text{ M.S.D}}{\text{Total V.S.D}} & 1 \text{ M.S.D} &= 1 \text{ cm} \\ & & 1 \text{ M.S.D} &= \frac{1}{10} \text{ cm} \\ &= \frac{0.1}{10} = \frac{0.1}{10} & &= 0.01 \text{ cm} \end{aligned}$$

$$= 0.01 \text{ cm}$$

Determination of zero error:

The type of instrument error found with the slide calliper is called zero error. These below are the three types of zero error.

- i) No zero error - If the zero of the main scale coincides with the zero of the vernier scale it is said as no zero error.
- ii) +ve zero error - If the zero of the vernier scale is on the right side of the main scale it is called as positive zero error.

- (iii) -ve zero error - If the zero of the vernier scale is on the left side of the main scale it is said to be negative zero error.

Procedure :-

- (i) The solid cylinder is to be presented in the form of diameter in such a way that both the jaws of the vernier callipers holds the solid cylinder.
- (ii) The solid cylinder placed in between the two jaws and it is changed a various no. of times to get the diameter.
- (iii) The readings are to be taken while doing each of the observation.
- (iv) The length of the solid cylinder is taken when the cylinder is placed in parallel to the slide calliper.
- (v) The procedure is repeated to get the length of the solid cylinder.
- (vi) The experiments are repeated and the observations are considered.

Observation :-

Tabulation-1
(External Diameter of the solid cylinder)

Sl No.	Least count in cm (L.C.)	Main scale reading in cm (M.S.R)	vernier coincidence (V.C.)	vernier scale reading in cm (V.S.R)	Total diameter (M.S.R + V.S.R)	Mean error in cm	zero error in cm	Corrected diameter in cm
1	0.01	1.2	9	0.09	1.29			
2	0.01	1.2	8	0.08	1.28			
3	0.01	1.2	8	0.08	1.28			
4	0.01	1.2	9	0.09	1.29			
5	0.01	1.2	9	0.09	1.29	0	0	1.288
6	0.01	1.2	9	0.09	1.29			
7	0.01	1.2	9	0.09	1.29			
8	0.01	1.2	9	0.09	1.29			
9	0.01	1.2	9	0.09	1.29			
10	0.01	1.2	9	0.09	1.29			

Tabulation - II (Length of solid cylinder)

SL. No.	Least Count in cms (L.C)	Main Scale reading in cms (M.S.R)	Vernier coincidence (V.C)	Vernier Scale reading in cms (V.S.R)	Observed Length in cms (M.S.R + V.S.R)	Mean error in cm	Zero Error in cm	Corrected Length in cms (D)
1	0.01	10.1	6	0.06	10.16			
2	0.01	10.1	4	0.04	10.14			
3	0.01	10.1	3	0.03	10.13			
4	0.01	10.1	4	0.04	10.14	0	0	10.145
5	0.01	10.1	4	0.04	10.14			
6	0.01	10.1	4	0.04	10.14			
7	0.01	10.1	5	0.05	10.15			
8	0.01	10.1	6	0.06	10.16			
9	0.01	10.1	4	0.04	10.14			
10	0.01	10.1	5	0.05	10.15			

Calculation:-

The diameter of the solid cylinder = 1.288

So, radius = $\frac{D}{2} = \frac{1.288}{2} = 0.644$

And, length of the solid cylinder (l) = 10.145

Hence,

$$\text{Volume (V)} = \pi r^2 l \text{ cm}^3$$

$$= \frac{22}{7} \times (0.644)^2 \times 10.145$$

$$= 3.14 \times 0.414736 \times 10.145$$

$$= 13.21153$$

Expt. No. _____

Date _____

Page No. _____

Result ~

The volume of the solid cylinder is found ~~1321153~~

Submitted by :-

NAME - Hitesh kumar Garraiyak

R. N - F20159001007


Class - Diploma 1st year

Group -

Branch - Civil

Date - 05-03-2021

Teacher's Signature _____


05/03/21

AIM OF THE EXPERIMENT :

To find the volume of a hollow cylinder using a vernier caliper.

APPARATUS REQUIRED :

Slide caliper

Hollow cylinder meter scale

THEORY SCALE :

Depth of the cylinder = Main scale reading + vernier scale reading = M.S.R + (V.C + L.C)

CALCULATION OF L.C :

(a) Standardization of main scale : $10 \text{ MSD} = 1 \text{ CM}$
 $1 \text{ MSD} = 0.1 \text{ CM}$

(b) Least count = $1 \text{ MSD} = 1 \text{ VSD}$
 (or) $10 \text{ VSD} = 1 \text{ MSD}$
 $10 \text{ VSD} = 0.9 \text{ CM}$
 So, $1 \text{ VSD} = 0.1 \text{ CM}$

Hence L.C = $0.1 \text{ cm} = 0.09 \text{ cm}$
 $= 0.01 \text{ cm}$

DETERMINATION OF ZERO ERROR :

The type of instrument error bound with the slide calliper is called zero error. These are the three types of zero error written below.

(i) No zero error :- If the zero of the main scale coincides with the zero of the vernier of is said as so.

Zero error

(ii) +ve zero error :- If the zero of the vernier scale is on the right side of the main scale it is said as the zero error.

(iii) -ve zero error :- If the zero of the vernier scale is on the left side of the main scale it is said as -ve zero error.

PROCEDURE :-

1. The L.C of the instrument is first calculated by taking the difference in the volume of 1 MSD and 1 VSD.
2. The zero error is found and by taking into the zero of the main scale and vernier scale when the jaws are closed.
3. The given hollow cylinder is introduced in the lower jaw and the external is found out.
4. The above procedure is repeated for at least 10 times to find the internal diameter in hollow cylinder.
5. The hollow cylinder is placed on the upper jaw to find the internal diameter of the given cylinder.

6. The procedure is repeated a least 10 times to get the internal diameter
7. The depth of hollow cylinder can be measured by the strips of the vernier caliper
8. The strips is introduced inside the hollow cylinder in such a way that the edge of the hollow cylinder.
9. The observation is taken at least 10 times.
10. All the readings for external diameter and depth of the cylinder are recorded in table.

THEORY

Let $D =$ External diameter of cylinder
 $d =$ Internal diameter
 $h =$ depth of cylinder

CALCULATION :-

$$\text{External diameter (D)} = 2.54 \text{ cm}$$

$$\text{Internal diameter (d)} = 2.23 \text{ cm}$$

$$\text{Depth (h)} = 7.52 \text{ cm}$$

Volume of hollow cylinder

$$= \frac{\pi}{4} (D^2 - d^2) h$$

$$= \frac{3.14}{4} \{ (2.54)^2 - (2.23)^2 \} 7.52$$

$$= 8.72 \text{ cm}^3$$

RESULT :- The volume of hollow cylinder is 8.72 cm³.

Submitted by Bhubendra Tudu


Roll No - F22159010014

Branch - Mining Engg

Group - A1, 2nd sem

Date - 09.05.23

Teacher's Signature


09/05/2023

AIM OF THE EXPERIMENT :-

Measurement of radius curvatures of convex surfaces of a watch glass.

APPARATUS REQUIRED :-

Spherometer

Watch glass

Base plate

THEORY :-

If R_1 is the radius of curvatures of convex surface of the given watch glass then

$$R_1 = \frac{d^2}{8h_1} + \frac{h_1}{2}$$

where h_1 = height of the convex surface from the base plate
 d = distance between any pair of fixed legs.

Determination of pitch and L.C of spherometer

10 division of linear scale = 1 cm

1 division of linear scale = 0.1 cm

1 complete rotation = 1 division of linear scale

Hence, pitch = 0.1 cm


$$L.C = \frac{\text{Pitch}}{\text{Total no. of division on the circular scale}}$$

$$= \frac{0.1 \text{ cm}}{100} = 0.001 \text{ cm}$$


$$\text{Height} = \text{Pitch Scale reading} + \text{circular scale reading} \\ = p$$

PROCEDURE :-


- (1) The pitch and L.C of the spherometer are determined.
- (2) The watch glass is placed with convex surface up on the base plate and over the convex surface the spherometer is placed.
- (3) The central leg is moved downwards gradually in such a way that it just touches the top protruding part of convex surface at this position the tip of central leg and its reflecting image just coincides.
- (4) The initial circular scale coincidence is noted at such position. The watch glass is then taken and the base and the spherometer is placed on it.
- (5) The central leg is gradually lowered down till it just touches the base plate. The no. of complete rotations is noted and then the circular scale is bound out when central leg touches the base plate.
- (6) All these readings are recorded on the table and the height of convex surface is calculated.
- (7) The above procedure is repeated for at least ten more observations and mean value of h is calculated.
- (8) To find 'd' the three legs of the spherometer are pressed on a graph paper at three different locations. The three points obtained from pressing are going to form an equilateral triangle. Three such figures are obtained from which mean value of 'd' is calculated.
- (9) The value of the radius of curvature is calculated by using formula given in the theory.



D_1



D_2



D_3

CALCULATION:-

Height of the convex surface = 0.2768 cm

$$d = 4.06 \text{ cm}$$

$$R_1 = \frac{d^2}{6h_1} + \frac{h_1}{2}$$

$$= \frac{(4.06)^2}{6 \times 0.27} + \frac{0.27}{2}$$

$$= \frac{16.48}{1.62} + \frac{0.27}{2}$$

$$= 10.367 \text{ cm}$$

RESULT:- The radius of curvature of convex surface is 10.367 cm

Submitted by Bhuberdra Tudu


F22157010019

Mining Engg, 5th sem

Group - A1

Date = 16.05.23

Teacher's Signature



AIM OF THE EXPERIMENT :-

To determine the radius of curvature of concave surface by using spherometer.

APPARATUS REQUIRED :-

- i) Spherometer
- ii) watch glass
- iii) base plate

THEORY

R_g is the radius of curvature of concave and concave surface of the given watch glass then

$$R_g = \frac{d^2}{6h_g} + \frac{h_g}{2}$$

h_g = Height of the concave surface from the base plate.

Determination of pitch and L.C of the spherometer

~~10 division of linear scale = 1 cm~~

~~1 division of linear scale = 0.1 cm~~

~~1 complete rotation = 1 division of linear scale~~

Hence, pitch = 0.1 cm

$$L.C = \frac{\text{Pitch}}{\text{Total no. division on the circular scale}}$$

$$= \frac{0.1 \text{ cm}}{100} = 0.001 \text{ cm}$$

Teacher's Signature _____

PROCEDURE :-

- i) The pitch and L_c of the spherometer are determined.
- ii) The spherometer is placed over the concave surface of the base plate.
- iii) The central leg is moved downwards gradually in such a way that it just touches the lowest part of the concave surface.
- iv) The initial circular scale coincidence is noted at such position.
- v) The central leg is gradually rotated upward till it touches and remains balanced on the base plate then the final circular scale is noted.
- vi) The above procedures are repeated for at least ten more observations and the mean value of h_s is calculated.
- vii) To find the value of d , the three legs of the spherometer is passed on a paper at three different positions. The three points obtained from pressing is joined to get an equilateral triangle. Three such triangles are obtained from which mean value of d is calculated.
- viii) The value of R_s is calculated by using the formula given in the theory.

Teacher's Signature _____

OBSERVATION:-

Table-1

Sl No.	Pitch in cm	L.C in cm	Initial scale reading	No. of rotation (N)	Final scale reading (F)	Difference 1~F	pitch scale reading	Circular scale reading L.C.R.D	Total height PSR + CSR	mean
1	0.1	0.001	22	6	1	31	0.6	0.031	0.631	
2	0.1	0.001	45	6	2	42	0.6	0.042	0.642	
3	0.1	0.001	35	6	3	32	0.6	0.032	0.632	
4	0.1	0.001	30	6	5	25	0.6	0.025	0.625	$\frac{6.935}{10}$
5	0.1	0.001	40	6	4	36	0.6	0.036	0.636	$= 0.6225$
6	0.1	0.001	37	6	10	27	0.6	0.027	0.627	
7	0.1	0.001	27	6	15	12	0.6	0.012	0.612	
8	0.1	0.001	20	6	20	2	0.6	0.002	0.602	
9	0.1	0.001	20	6	9	11	0.6	0.011	0.611	
10	0.1	0.001	25	6	10	15	0.6	0.015	0.615	

Table-2

Sl No	D ₁ in cm	D ₂ in cm	D ₃ in cm	$D = \frac{D_1 + D_2 + D_3}{3}$	mean of
1	4	4.5	4.2	4.23	$\frac{12.69}{3}$
2	4.2	4	4.5	4.23	$= 4.28$
3	4.5	4.2	4	4.23	

Teacher's Signature _____

CALCULATION:-

Height of the concave surface = 0.6225 cm
 $d = 4.23$ cm

$$\begin{aligned}
 R_2 &= \frac{d^2}{6h_2} + \frac{h_2}{2} \\
 &= \frac{(4.23)^2}{6 \times 0.62} + \frac{0.62}{2} \\
 &= \frac{17.89}{3.72} + 0.31 \\
 &= 4.80 + 0.31 \\
 &= 5.11 \text{ cm}
 \end{aligned}$$

RESULT :- The radius curvatures of concave surface is 5.11 cm

Submitted by Bhubendra Tudu


F22159010014




mining engg, 8th sem

Group - A1

Date - 30.05.23

Teacher's Signature





AIM OF THE EXPERIMENT :-

To study variation of time period of a simple pendulum of a given length by taking bobs of same size but different masses and interpret the result.

APPARATUS REQUIRED :-

A spherical metallic bob with a bob, 2 meters long thread, spirit level, clamp stand, slide callipers, metre, scale, stop watch.

THEORY :-

The time period of a simple pendulum is given by

$$T = 2\pi\sqrt{\frac{l}{g}} \quad \text{or} \quad T^2 = 4\pi^2 \frac{l}{g}$$

$$\Rightarrow \boxed{g = 4\pi^2 \frac{l}{T^2}}$$

where, l = effective length of the simple pendulum
 g = acceleration due to gravity at the place of observation.

PROCEDURE :-

- (1) The value of (2π) is determined by subtracting the value of radius of the bob from $2\pi r$. The values of r and $2\pi r$ are calculated with the help of slide calliper.
- (2) A thread is taken and small ink marks are given at places where the length of the thread when added to the value of r gives effective length a whole number like 40 cm, 50 cm, 60 cm respectively.

- (3) The thread is now introduced between the split- cork and clamped in such a way that the point of suspension (lower surface of the cork) from the centre of the bob is 120 cm.
- (4) A chalk mark is given on the edge of the table as mean position of pendulum. Two chalk marks are drawn about 3 cm on either side of the mean like to a part limit the amplitude of oscillation to 4° .
- (5) The pendulum is now allowed to oscillate about the mean position in a vertical plane without spinning. For this the bob is drawn parallel to the edge of the tube through a small distance and allowed to go for oscillation.
- (6) The time period for 20 oscillation is noted by stop watch this process is repeated thrice and the mean value of the period is calculated - The time period 'T' of simple pendulum for each bob for effective length 120 cm.
- (7) The pendulum is suspended from one effective length and the mean time for 20 oscillations is calculated. The time period 'T' of simple pendulum for each bob for same effective length obtained by dividing the mean value by 20.
- (8) Different bob of same size but different masses are taken and similarly time period T is found out.
- (9) Graph are drawn taking effective length along X-axis and T^2 along y axis straight line passing through the origin are obtained.
- (10) From graph $2/\pi^2$ is calculated in each case and acceleration due to gravity (g) is found out.

RESULT :- Hence, the acceleration due to gravity 'g' of simple pendulum is 973.95 cm/s^2 .

Submitted by Bhubendra Tudu

F22159010014

Mining Engg, 2nd sem

Group - A₁

Date - 31.05.23

Teacher's Signature _____


21/05/23

AIM OF THE EXPERIMENT :-

To find the angle of prism

Apparatus Required :-

- i) Drawing Board
- ii) Fixing pins
- iii) Prism
- iv) protractor
- v) Drawing pins

THEORY

The angle between two surfaces of a prism is known as refracting angle or angle of prism. The angle of prism is 60° i.e. $\angle A = 60^\circ$.

PROCEDURE :-

1. Fix the sheet of paper on the drawing board by stretching its edges.
2. Place the prism on the paper and draw the outline of it.
3. Draw a parallel line on both the sides of prism from the mid points of the prism.
4. Place two pins on either side of the prism and by looking through the prism the other two pins are fixed which is given the refraction. It looks like in straight line.
5. Repeat the process for the other side of prism.
6. The experiment is to be repeated for at least 5 times and the mean should be calculated.

Expt. No. 08

Date 06.06.93

Page No. 31

OBSERVATION :-

Sl. No.	Angle of prism L.A	Mean
---------	-----------------------	------

RESULT :-

Therefore the Angle of prism to be 60.9° .

Name - Bhubendra Tripathi

Roll no - FG2159010014

Branch - Mining

Semester - 2nd

Date - 06.06.23

AIM OF THE EXPERIMENT :-

To determine the angle of minimum deviation by I-D

APPARATUS REQUIRED :-

- i) Prism
- ii) Drawing pine
- iii) Fixing pins
- iv) protractor
- v) Drawing board

THEORY :-

A prism is an optical element. Its polished flat surfaces that reflect light.

- when light travels from one medium to another medium it is reflected at the second medium at an different angle.
- The degree of bending of the light path depends on the angle that the incident beam of light makes with the surface of the prism. μ is the ratio between the refractive index of the two medium, this is called Snell's law i.e. $\mu = \frac{\sin i}{\sin r}$

Where (μ) is the refractive index of the material of the prism (i) is the angle of incidence (r) is the angle of refraction.

→ The angle of minimum deviation can be obtained from the graph of (α) be the angle of minimum deviation then the refractive index (μ). At the material of the prism is calculated using the formula.

$$\mu = \frac{\sin(A + \frac{D_m}{2})}{\sin \frac{A}{2}}$$

PROCEDURE :-

- i) A paper is fixed in the drawing board placed on the table.
- ii) place the given glass prism on the paper, using the pencil work the outline A, B, C of the prism in the paper.
- iii) Remove the prism using the scale pencil hermet 'No' is drawn, using the protractor measure an angle 30° the normal.
- iv) Another line (incident) is drawn the angle 30° with the normal 'No'.
- v) Two pins P_1 and P_2 are fixed on the line, the prism is reflected on the outline. By looking the pins from AC of the prism then other pins P_3 and P_4 fixed that P_1, P_2, P_3, P_4 are in one line. Remove the pins how extend the incident line in one line forward direction and the refracted line in the backward direction, they meet at a point let the point be P.
- vi) Using the protractor measure the angle at the point this is the angle at the point P this is the angle of deviation repeat the experiment for different value.
- vii) Repeat the experiment for different angle of incidence $i = 30^\circ, 35^\circ, 40^\circ, 45^\circ, 50^\circ, 55^\circ$ of the corresponding angle of deviation measured.

- viii) Draw a graph angle of incidence (i) an angle of deviation along with the 'y' axis.
- (ix) The angle of deviation corresponding to the lowest of the curve is the angle of minimum deviation.
- (x) Angle of the prism 'A' is directly measured from the outline of the prism using the protractor calculate the refractive unit of the material of the prism using the formula.

OBSERVATION:

Sl. No	Angle of incidence	Angle of deviation	Angle of minimum deviation
1.			
2.			
3.			
4.			
5.			
6.			

$$\mu = \frac{\sin \left(A + \frac{D_m}{2} \right)}{\sin \frac{A}{2}}$$

Date 23.06.23

Page No. 35

Expt. No. 09

RESULT :-

The angle of minimum deviation by I-D is 1.92

Submitted by Bhubendra Tudu

Roll No. - F22159010014

Branch - Mining Engrg

Semester - 2nd Sem.

Group - A1

Date - 23.06.23

Teacher's Signature


23/06/23

AIM OF THE EXPERIMENT:

To trace lines of force due to a bar magnet with north pole pointing south and locate the neutral points.

APPARATUS REQUIRED:-

- i) Bar magnet
- ii) Compass box
- iii) Drawing board
- iv) Scale
- v) pencil
- vi) fixing pins.

THEORY:-

A magnetic line of force is a closed curve that indicates the path along which a unit north pole can move to its force to do so within a magnetic field.

PROPERTIES:-

- (i) They always start from the north pole and end on south pole and continues throughout the body of the magnet.
- (ii) The tangent at any point on the line of force gives the direction of magnetic field at that point.
- (iii) The no. of lines of force per unit area being perpendicular at P is directly proportional to the magnitude of the strength of the field.


- iv) The two lines of force do not intersect with each other. As it shows at the point of intersection it will give tangents which will show two divergent sections which is not possible.
- v) Neutral points are the points in the region around a magnet where its magnetic field is neutralised by the earth's magnetic field i.e. the resultant force is at this point.

Conclusion:-

In the above manner the line of force due to a bar magnet with north pole pointing south and the neutral points are located.

submitted by Bhubendra Tudu
P22159010014
Mining Engg, 2nd sem
A1 group
Date - 23.06.23

Teacher's Signature


23/06/23

AIM OF THE EXPERIMENT:

To trace the lines of force due to a bar magnet with north pole pointing north and locate neutral points.

APPARATUS REQUIRED:—

- i) Bar magnet
- ii) Compass box
- iii) Drawing board
- iv) Scale
- v) pencil
- vi) Fixing pins

THEORY:—

Magnetic lines of force:

A magnetic line of force is a closed curve that indicates the path along which a unit north pole free to do so within a magnetic field.

PROPERTIES:

- 1) They always start from the north pole and end on south pole and continuous throughout the body of the magnet.
- 2) The tangent at any points on the lines of force gives the direction of magnetic field at that point.
- 3) The no. of lines of force per unit area (area being perpendicular) is \propto directly proportional to the magnitude of the strength of the field.

4) The two lines of force do not intersect with each other. If it shows at the points of intersection it will give tangents which will show two directions which is not possible.

5) Neutral points in the region around a magnet when its magnetic field is neutralised by the earth's magnetic field i.e. the resultant force is zero at these points.


Conclusion:

This is the line of force due to a bar magnet with north pole pointing north and the neutral points are located.

Done

submitted by Bhubendra Tudra
F22159010014
Mining Engg/
2nd sem
A1 group
Date - 27.06.23

Teacher's Signature


27/06/23